

- Buffer solutions
- Cleaning agent
- Probe stand
- Probe extension cord
- Dust caps
- Power supply
- Batteries

Contact:

Sentron Europe BV
Aan de vaart 3
9301 ZH Roden
The Netherlands

info@sentron.nl
tel +31 50 501 38 00
fax +31 50 501 68 34

| Item | Description | Order code | Details |
|------|-------------|------------|---------|
|------|-------------|------------|---------|

Buffer solutions



Sentron buffer solutions are N.I.S.T. based, high quality reference fluids.

The twin neck bottle provides a convenient and easy to use way of providing the right amount of buffer without waste or chance on contamination.

| | |
|----------|---|
| A116-002 | pH 2.00 buffer 0.5l +/- 0.02 pH N.I.S.T @ 25°C |
| A116-004 | pH 4.00 buffer 0.5l +/- 0.01 pH N.I.S.T @ 25°C |
| A116-007 | pH 7.00 buffer 0.5l +/- 0.01 pH N.I.S.T @ 25°C |
| A116-010 | pH 10.00 buffer 0.5l +/- 0.01 pH N.I.S.T @ 25°C |
| A116-012 | pH 12.00 buffer 0.5l +/- 0.05 pH N.I.S.T @ 25°C |

Cleaning agent



Cleaning agent for removal of enzymes and proteins.

Alconox Terg-A-Zyme is a concentrated powder and can be prepared by adding (deionized) water. One package makes approx. 200 liters of detergent.

| | |
|----------|--------------------------|
| A118-011 | 1.8kg of cleaning agent. |
|----------|--------------------------|

Probe stand



The probe stand can fit probes and includes 2 adapter rings for mounting probes with various shapes and dimensions.

The probe stand has a solid weighted base and an easy adjustable arm to which the probe cable(s) can be attached.

| | |
|----------|------------------------------------|
| A100-001 | Probe stand incl. 5 adapter rings. |
|----------|------------------------------------|

Probe extension cord



Probe extension cord of 1.5 m.

The extension cord provides additional length for those situations where extra distance between probe and meter is required.

| | |
|----------|-----------------------------|
| A109-010 | Probe extension cord 1.5 m. |
|----------|-----------------------------|

| Item | Description | Order code | Details |
|---|--|------------|---|
| <p>Dust caps</p>  | <p>A set of dust caps for both meter and probe connector.</p> <p>The dust caps completely seal and cover the connectors against dust and moisture when the probe and/or meter is not in use.</p> | A102-002 | <p>A set of two dust caps for both the meter and the probe connector.</p> |
| <p>Power supply</p>  | <p>Power supply including five world plug adapters and a USB power cord.</p> <p>Adapters included are: UK, US, EU, Australia, Middle- and Far East.</p> | A102-006 | <p>Power adapter set <i>in:</i> 100-240 VAC 50/60Hz, 0.2A <i>Out:</i> 5VDC , 1A</p> |
| <p>Batteries</p>  | <p>Replacement battery set of six AA 1.2V batteries for SI meters.</p> <p>High quality - 2400mAh NiMH - rechargeable batteries which are recommended for use in SI meters.</p> | A102-011 | <p>Battery set of six batteries AA 1.2V, 2400 mAh NiMH (ANSI-1.2H2 or IEC-HR6)</p> |